# DRY FLUOROIMMUNOASSAY DEVICE IVD DEVICE

MEDICINES ONLINE

#### RRP £1,998

#### **Key benefits**

- Touch screen.
- Flu or COVID result.
- Up to 69 different tests.
- Purchase or rentals available.
- Intelligent QR code management.
- LIMS compatible data connectivity.
- Results in 15-30 minutes no need to send to a lab.
- Ultra-high Sensitivity and excellent precision.

### **Overview**

• The instrument represents advanced technologies of the same type and can test relevant parameters of the human body conveniently, accurately and quickly when used cooperatively with fluoroimmunoassay reagent

strips.

- KD-120P Dry Fluoroimmunoassay Analyzer uses a LED as the excitation light source. The light from the LED hits a test device that has been inserted into the meter. This causes the fluorescence dye in the test device to give off energy. The more energy the fluorescent dye gives off, the stronger the signal.
- After a buffer-mixed sample is added to the test device, the test device is inserted into the KD-120P Dry Fluoroimmunoassay Analyzer.

## **Specifications**

ID CHIP

 $(\epsilon)$ 

Item	Parameter	
Excitation light source	LED	
Wavelength	365nm	
Sample type	Whole blood, serum, plasma, urine and faeces (subject to specification of immunofluorescence assay ki	
Test speed	240 times/hour	
Power source	A.C.100V-240V 50Hz/60Hz Host: DC12V 5A	
Input power	60 VA	
Background number	Fluorescence value should be $\leq 5$	
Repeatability	The variable coefficient of fluorescence value of fluorescence test cards of different intensities repeatedly measured by the instrument CV≤2.0%	
Stability	The relative deviation between the 4th hour and the 8th hour after the instrument is working stable, and the test result when the device just works stably is limited to $\pm 3.0\%$ .	
Accuracy	When the instrument tests the calibrated fluorescence test card, the deviation between the measured fluorescence value and the calibration value of the fluorescence test card is ±3.0%.	
Measuring time	The time from inserting the fluorescence test card to displaying the test result should not exceed 15s.	
Linearity range	The linearly dependent coefficient of the fluorescence value and concentration read by the instrument $r \ge 0.990$	

### Test List

	TSH
	TT3
Thyroid Hormone	TT4
	FT3
	FT4
Renal injury	CysC
	NGAL
	MAU
	RBP-4
	β-MG
	cTnl
Cardiac	NT-proBNP
Cardiac	D-Dimer
	Lp-PLA2
	CRP/hs-CRP
	SAA
Inflammation	CRP+SAA
	IL-6
Diabetes	Insulin
Diabetes	C-peptide

	β-HCG
	АМН
	LH
	SHBG
Hormones	FSH
	PRL
	Progesterone
-	Testosterone
	Oestradiol
	Vitamin D
Nutritional metabolism	Vitamin B12
	FERR
	AFP
	CEA
	NSE
Tumour markers	CA-125
	CA15-3
	β-HCG
	HE4
	CA19-9

\*Some tests are part of Medicines Online research program. Please provide feedback.



## MEDICINES ONLINE<sup>®</sup>

27 Old Gloucester Street | London | WC1N 3AX