



INSTRUCTIONS URINE STI HOME TEST KIT

PLEASE READ ALL INSTRUCTIONS BEFORE COLLECTING YOUR SAMPLE

YOUR KIT CONTAINS



Sample container



Return bag

BEFORE TAKING YOUR SAMPLE

ORDER NO.	DAY	CUSTOMER SERVICE	REC.	LAB STAGE	DEL. SERV.	RESULT
12345	01/12/2023	01/12/2023	01/12/2023	01/12/2023	01/12/2023	01/12/2023

Login to your account on www.medicinesonline.org.uk and activate your sample by clicking on **SAMPLE ACTIVATION AND RESULT**, then select your test day, enter **Royal Mail Tracking Number**, and **tube number (number under the barcode)** then click on save.

TAKING YOUR TEST

1.



Wash your hands and open the sterile collection sample.

2.


Urinate a small amount into the toilet. Place the sample container under stream and continue to urinate and collect specimen. About half the sample container is enough.

3.



Screw the lid on the sample container tightly. Do not touch the inside of the sample container or lid.

4.



Place the sample container into the return bag. And post the bag in a Royal Mail Priority Post box. To find the nearest Priority Post Box please download the Royal Mail app or go to www.royalmail.com/services-near-you/